wau

LED HAND MASK

PHOTOTHERAPY HAND MASK

USER GUIDE

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Table of Contents

Technical Specifications	26
Purpose	27
Device components	27
Operating Procedures.....	28
Usage Method	30
Precautions	31
Cleaning	32
Storage.....	33
Charging	33
Disposal.....	34
Warranty.....	34


Technical Specifications

Parameter	Value
Model	FM-05
Materials	ABS plastic, silicone
Dimensions	385x185x52 mm
Weight	326 g
Power source	Li-ion battery
Battery capacity	1500 mAh
Operating voltage	3,7 V
Supply voltage	5 V
Modes	Automatic and manual control
Light types and wavelengths	Red (630 nm) Blue (415 nm) NIR (Near Infrared light) (830 nm) Mixed (red + blue)
Smart-timer	Automatic control mode — 15 minutes Manual control mode — 10 minutes
Charging time	1,5 hours
Max continuous use time	40 minutes
Irridation intensity	Red — 4,9 mW/cm ² Blue — 28 mW/cm ² NIR (Near Infrared light) — 10,1 mW/cm ² Mixed (red + blue) — 15 mW/cm ²
Number of LEDs	106 lamps

Purpose

LED HAND MASK is a hand mask designed for deep restoration and rejuvenation of the skin of the hands, wrists, and forearms. The red spectrum is aimed at activating collagen synthesis, which helps restore skin elasticity and eliminate age-related pigmentation. The blue and purple light modes aim to reduce inflammation and provide an antiseptic effect, while Near Infrared (NIR) technology is designed for a deeper impact on the tissues. The use of safe food-grade silicone and a wireless controller makes the professional care procedure as comfortable and accessible as possible at home.


Device components

1. LED HAND MASK
2. Remote control
3. Type C-Type C charging cable
4. Fastening straps

Remote control

1. Power button / light selection
2. Auto and manual mode switch button
3. Automatic control mode
4. Manual control mode
5. Red light mode
6. Blue light mode
7. Near Infrared (NIR) light mode
8. Mixed light (red + blue) mode


Package contents


Mask x1
Fastening straps x2
Type C-Type C charging cable x1
User manual x1, storage bag x1

Operating Procedures

1. Attach the fastening straps to the mask. Place the mask face up. Thread the right strap (its front side is completely covered with velcro, and the back with silicone) into the right hole on the mask. Fasten one end of the strap to its body.

2. Repeat these steps with the left strap. Its front side has velcro, but is mostly covered with silicone. Thread the left strap into the hole on the left side of the mask. Fasten one end of the strap to the velcro.

3. Put the mask on your hand. Fasten the left strap to the right strap.
4. Short press the Power button. The mask automatically turns on in manual control mode with red light.
5. Using the mode switch button, select your desired control mode: automatic (lights change automatically) or manual mode (1 light stays on throughout the procedure).
6. After selecting one of the control types, an LED indicator will light up next to the mode symbol A (auto — automatic) or M (manual — manual) on the remote control.
7. In manual control mode, use the Power button / light selection to select the desired light for the phototherapy session.
8. The red light automatically starts first (symbol on the remote control — 1). Then, by pressing the button, the light switches to blue (symbol on the remote control — 2). Pressing again starts the Near Infrared light (symbol on the remote control — 3). Another press of the button turns on the mixed light (symbol on the remote control — 4).
9. In manual control mode, a 10-minute smart-timer operates; after it expires, the mask will turn off automatically.
10. When switching modes, the smart-timer will restart.
11. In automatic mode, the lights will change automatically: 5 minutes — red light, 4 minutes — blue light, 3 minutes — Near Infrared light, 3 minutes — mixed light.
12. The total session time in automatic mode is 15 minutes. Afterwards, the mask will turn off automatically.
13. To turn off the mask, press and hold the Power button.

Usage Method

Sensitivity test

[Important] You must carry out a sensitivity test before your first treatment.

- Place the device on the inside of your forearm.
- Turn on the mask and run at least a 10-minute phototherapy session.
- If you feel heat or discomfort, stop using the device immediately.
- After completing the test, wait at least 6 hours. If the skin on the inside of your arm is red and/or itchy, the LED HAND MASK is not suitable for you. Contact the organization authorized to handle consumer claims.

After completing the sensitivity test:

- Cleanse the skin of your hand before starting the session, do not apply skincare products.
- Attach the fastening straps to the mask.
- Put the mask on your hand, fasten the straps.
- Select your desired control mode: automatic or manual.


The following light types are available in LED HAND MASK:

Red (630 nm)
Stimulates blood circulation, reduces inflammation, heals microcracks, and lightens pigmentation spots.

Blue (415 nm)
Normalizes skin oiliness, destroys bacteria, fights breakouts, and soothes.

Near Infrared (NIR) (830 nm)
Penetrates into deep tissue layers, accelerates cell renewal, increases elasticity, and tightens the skin of the hands.

Mixed (red + blue) (630 + 415 nm)
Combines the power of two spectrums for a powerful anti-aging, antibacterial, and restorative effect.

- Run a phototherapy session. In manual control mode, it is not recommended to exceed the 10-minute usage limit.
- After the phototherapy session, it is recommended to apply skincare and moisturizing products.

Frequency and duration of use:

The mask is equipped with a smart-timer, in automatic control mode the session will end automatically after 15 minutes, in manual control mode the session will automatically end after 10 minutes. It is not recommended to use the device longer than the specified time.

If you are using the phototherapy mask for the first time, we recommend following the schedule below.

Week	Frequency	Manual control	Automatic control
Week 1	Every 3 days	10 minutes	15 minutes
Week 2	Every 2 days	10 minutes	15 minutes
Week 3	Every 2 days	10 minutes	15 minutes
Week 4	Every day	10 minutes	15 minutes

To achieve a visible effect, it is recommended to use the LED HAND MASK for at least 4 weeks.

Precautions

- Before the first use, carry out a sensitivity test to light (see the «Usage Method» section).
- The mask may only be used for its intended purpose in accordance with the user manual.
- Do not use the mask on skin areas with open wounds, ulcers, deep cracks, fresh postoperative scars, or signs of severe inflammation.
- Consult a physician before using the mask during pregnancy and lactation.
- Do not look directly at the lit diodes while the device is operating.
- Do not use or store the mask near heaters, open flames, under direct sunlight, or in a car in hot weather.
- Do not bend the mask too much, do not squeeze it, and do not place heavy objects on it, as there are flexible microcircuits and LEDs inside the mask — this can cause mechanical damage.
- Do not attempt to disassemble or modify the device yourself. This can lead to damage, electric shock, and will void the warranty.

Warnings:

- It is prohibited to use the device if you suffer from diseases causing increased photosensitivity (e. g., lupus, photodermatitis), or take medications/cosmetics that increase photosensitivity (retinoids, certain antibiotics, citrus essential oils).
- Stop using the device if you experience discomfort or worsening of your condition.
- Do not use the mask while it is plugged in and charging. Always disconnect the power cable before the procedure.
- Do not use the mask if it is damaged, contact the organization authorized by the manufacturer to handle consumer claims.
- The device is not waterproof. Do not use it in the bathroom, near sinks or other water sources. Do not touch the mask, charging cable, or adapter with wet hands.
- Clean the mask only after turning it off and unplugging it. Do not immerse the device in water and do not wash it under running water.

Cleaning

1. Clean the inner surface of the mask after each use.
2. Turn off the mask before cleaning.
3. Do not use high-temperature or aggressive disinfectants.
4. Do not use benzene, thinners, butyl acetate, or other flammable substances for cleaning.
5. Do not spray water directly onto the surface of the mask.

Surface cleaning:

1. Wipe both sides of the mask with a soft, lint-free cloth dampened with a mild water-based cleaning solution.
2. Place the mask on a clean, dry surface until it is completely air-dried.

Visual inspection:

1. Inspect the device surface for any visible dirt. If contamination remains, proceed to the next step.

Deep cleaning (for stubborn stains):

1. Carefully wipe the silicone surface twice with a 70%–80% ethanol solution (3 minutes each time), then wipe thoroughly dry with a damp or dry cloth.
2. Allow the mask to air-dry completely before use.

Storage

Store the LED HAND MASK:


1. Away from flammable and explosive materials (sprays, oils, gas, etc.).
2. At room temperature in a dry place.
3. Avoid storing in dusty areas, in direct sunlight, or in high-humidity conditions.

Charging

When the battery level is low, the light indicator under the number 1 on the controller will start to flash. After six flashes, accompanied by beeps, the device will turn off automatically.

During charging, use only the original Type-C charging cable included in the package. It is strictly prohibited to use the hand mask while it is plugged in and charging. Using a power adapter with a voltage below 5V (V) will cause a charging failure. Using an adapter with a voltage above 5V (V) may damage the remote control.

To charge, insert one end of the Type C-Type C charging cable into the charging port on the mask controller, and connect the other end to the adapter.

The charging process is accompanied by the sequential flashing of the LEDs, which reflect the current power level of the built-in battery:

0-25%: The first LED indicator flashes, the rest are off.

25-50%: The first LED stays lit continuously, the second indicator flashes.

50-75%: The first two LEDs stay lit continuously, the third indicator flashes.

75-100%: Three LEDs stay lit continuously, the fourth indicator flashes.

100% (Charging complete): All 4 LED indicators stay lit continuously.

Disposal

At the end of the device's service life, its components (especially the lithium-ion battery and electronic parts) may harm the environment if disposed of with regular household waste.

1. Do not dispose of the device in a regular waste bin along with food waste.
2. Dispose of the device in accordance with local legal regulations.
3. Take the device to specialized electronics and battery collection points for further sorting and recycling.

Warranty

The WAU brand LED HAND MASK phototherapy hand mask is covered by a 1-year warranty. During the warranty period, the product is subject to replacement or repair if technical defects are discovered. This warranty does not cover malfunctions caused by improper operation or transportation, nor does it apply to products after unauthorized repair.

LED HAND MASK

PHOTOTHERAPY HAND MASK

Model: **FM-05**

Materials: **ABS plastic, silicone**

Dimensions (H×W×D): **385x185x52 mm**

Weight: **326 g**

Voltage: **3,7 V**

Power: **5,5 W**

Origin: **PRC**

Authorized manufacturer representative: WAU Global Limited.
63A College road, Northcote, Auckland, 0627, New Zealand.
e-mail: info@wauglobal.com

Important Notice: Important Notice: To ensure continuous product improvement, the manufacturer reserves the right to modify design and technical specifications. The latest version of this manual is available at: wauglobal.com