wau

MEI1

FACE AND NECK EMS MASSAGER

USER GUIDE

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RoHS

Table of Contents

Technical Specifications	24
Purpose	24
Device Components.....	25
Operating Procedure	26
Usage Method	26
Cleaning & Maintenance.....	32
Storage.....	32
Charging	32
Precautions	33
Disposal.....	33
Warranty	34


Technical Specifications

Parameter	Value
Model	BPWE-881
Materials	ABS plastic
Dimensions (HxWxD)	140x80x10 mm
Weight	78 g
Power	4,4 W
Operating Voltage	5 V
Power Source	Built-in lithium battery
Battery Capacity	800 mAh
Technologies	EMS, Wave vibration, Ion import (iontophoresis), Heating (42 °C)
Current	1,4 mA
Vibration	8 500 rpm
Smart Timer	10 minutes
Charging Time	≤ 2 hours
Battery Life (Full Charge)	80 minutes
Waterproof Rating	IPX3

Purpose

The Face and neck EMS massager MEI1 by WAU combines traditional Gua Sha massage with 4 advanced technologies: EMS, vibration, heating, and iontophoresis. Featuring two movable conductive spheres, the MEI1 deeply targets muscles and lifts facial contours, while its side surfaces combat puffiness and promote lymphatic drainage. The thermal effect combined with iontophoresis helps skincare products penetrate deeper, significantly enhancing their effectiveness.

Device Components

1. Silicone protective cover

2. Roller massage head

3. Power / Mode button

4. Mode indicator light ring

5. USB port rubber plug

6. Scraper / Ion import panel

Package Contents

Massager x1

Storage pouch x1

Type-C charging cable x1


User guide x1

Operating Procedure

1. Long press the power button for 1.2 seconds to turn the device on or off. A power-on/off beep will sound.

2. First short press (Default Mode): Activates the EMS + Wave Vibration + Ion Import combination. In this mode, the indicator light ring around the button glows white. This mode is designed for basic lymphatic drainage, reducing morning puffiness, toning muscles, and enhancing the penetration of daily skincare products.

3. Second short press (Comprehensive Mode): Activates the EMS + Wave Vibration + Ion Import + Deep Heating combination. In this mode, the indicator light ring around the button switches to red. This mode is designed for intensive anti-aging lifting, deep skin nourishment, relaxing tensed muscles, and triggering cellular rejuvenation processes through warming thermal massage.

4. The device operates continuously for 10 minutes. Switching modes does not reset the smart timer. The device will turn off automatically after 10 minutes.

5. The device will automatically turn off to protect the circuit if no action is taken for 2 minutes.

Usage Method

Facial and Neck Care:

1. Cleanse your face as part of your daily skincare routine.
2. Generously apply a conductive gel, hydrating serum, or essence that is strictly water-based. To use this device safely and effectively, your skincare product must not contain oils or acids (AHA, BHA, PHA acids, retinol, or high concentrations of Vitamin C).
3. Slide the device smoothly over your skin at a steady pace.

Selecting the Working Surface:

Move the device smoothly, at a constant speed, and without applying excessive pressure. The 360° design allows all sides of the housing to be used as active working surfaces:

- **Roller Head (Top Part):** Provides a deep, kneading 3D massage. It lifts the skin fold, sculpts the facial contours and cheekbones, and reduces puffiness.
- **Side Edges (Scraper / Ion Import Panels):** Designed for anatomical Gua-Sha-massage. They relieve muscle tension and ensure firm skin contact for EMS and deep product absorption.
- **Bottom Rounded Protrusion:** Intended for targeted care on localized areas (between the brows, sides of the nose) and acupressure massage on trigger points around the collarbones.

Recommended Massage Techniques:

Targeted Care:

Use the bottom rounded protrusion. Apply gentle pressure to trigger points, moving from the center outward along the upper and lower lines of the collarbones. This helps relieve accumulated fatigue, tightness, and muscle tension.

Lifting Care:

Press the smooth side panel (scraper) firmly against the skin.

Neck Area:

Begin the procedure on the neck to stimulate lymphatic flow. Smoothly guide the scraper downward, from the chin toward the collarbones, strictly avoiding the thyroid gland area.

Facial Area:

Massage moving upward from the jawline to the hairline on the forehead. Work on each facial zone by gliding the massager smoothly from the center of the face outward (toward the temples and ears), working against gravity to achieve a lifting effect.

**Roller Care (3D Sculpting):**

Roll the textured rollers over the skin for a deep V-shaped lift, calming, and muscle relaxation. Position the contours of your face so that the skin fold is caught between the two rollers.

Neck Area:


Roll the rollers up and down along the side muscles of the neck to relieve spasms and improve lymphatic drainage.

Facial Area:

Perform the massage by moving from the chin upward to the forehead hairline. Treat each zone by moving from the center of the face outward. For example: place the rollers on the lower third of the face and roll them along the jawline from the center of the chin to the earlobes, shaping a sharp jawline and defining the cheekbones.

Body Care:

Move the device smoothly over the skin along the massage lines. Thanks to its ergonomic shape, the device perfectly adapts to the anatomical contours of the body, combining deep massage technologies to relieve muscle tension, eliminate stagnant lymph, and remove excess extracellular fluid.

**Neck and Shoulder Care:**

For maximum results, combine the device's working surfaces: use the rollers for intensive kneading of muscle knots, and the side edges (scraper) for lymphatic drainage. Perform the massage along the following lines:

Downward Movements:


Smoothly guide the device from the base of the skull down to the shoulder joint to drain excess fluid and relieve spasms.

Vertical Movements:

Glide the rollers strictly up and down in the sternum area (or along the spine, avoiding pressure on the bones) to eliminate muscle stiffness.


Zigzag Movements:

Move the device in a zigzag pattern from the center outward for deep relaxation of trigger points.

Arm Care:

This massage targets active lifting, tightening sagging areas, and sculpting arm contours. Make smooth vertical movements alternately downward — from the armpit to the inner elbow — and then back upward in the opposite direction.

Leg Care:

Sculpts the silhouette, increases skin elasticity, and relieves tension. Guide the device strictly upward in the direction of lymphatic flow (from the ankles to the knees and thighs):

Massage Rollers:

Deeply knead tissues, improve blood circulation, reduce volume, and relax muscles.

Side Edges (Scraper):

Provide lymphatic drainage, eliminate excess fluid, reduce puffiness, and relieve «heavy legs» syndrome.

Frequency and Duration of Use

The massager is equipped with a smart timer; the session will end automatically after 10 minutes. It is not recommended to use the device longer than the specified time.

If you are using the EMS massager for the first time, we recommend following the schedule below:

Week 1: every 3 days for 5–10 minutes. Start with a 5-minute session, gradually increasing its duration.

Week 2: every 3 days for 10 minutes.

Week 3 and beyond: every 2 days for 10 minutes.

To achieve a visible effect, it is recommended to use the MEI1 for at least 4 weeks.

Cleaning

1. Turn off the power and switch off the device before cleaning.
2. After each use, wipe the rollers and the body with a damp cloth and keep them dry.
3. The device is not waterproof; do not use it in water, do not wash it under running water, and protect it from moisture.

Storage

1. For long periods of non-use, store the device fully charged and recharge the battery every 3 months.
2. Store in a cool, dry place, away from direct sunlight and high temperatures.
3. Store the device in a dust-free environment with low vibration levels.
4. Keep away from extreme cold, extreme heat, and open flames.
5. Keep out of reach of children.

Charging

1. Use only the original or recommended power adapter for this device.
2. Place the device in a well-ventilated area during charging; do not cover it with any objects.
3. If you notice unusual heating, smoke, or other abnormal conditions during charging, disconnect the power immediately.
4. Do not use the device while it is charging.
5. Please charge the device before first use to ensure efficient performance.
6. If there is no response when turning on the device, the built-in lithium battery may be depleted; charge it before trying again.

Low Battery Indicator:

When the battery is close to full depletion, the device will alert you: the indicator light ring will flash 10 times, accompanied by 10 beeps. After this, the device must be placed on charge.

Charging Process:

While connected to the power supply and charging, the indicator light ring will pulse with a breathing red light.

Fully Charged:

Once the battery is 100% charged, the indicator light will stop flashing and remain solidly lit, signaling that the device is fully ready for use.

Precautions

- This device is intended for personal home use only and is not a medical device.
- This device is intended only for its direct purpose as specified in this manual.
- Pregnant women should use the device only after consulting a doctor.
- Use with caution during menstruation or lactation.
- Use with caution if you have telangiectasia (rosacea/couperose).
- Use with caution if you are undergoing dental treatment or have dental fillings.
- Use with caution if you have a metal allergy.
- Use with caution if you are undergoing medical treatment, running a fever, or have hypertension.
- Do not use on cuts, plastic surgery areas, purulent or inflamed skin, itchy or red areas, or on the eyeballs or eyelids.
- Do not use on metal implants or anatomically sensitive areas (eyeballs, inner ears, bones, joints).
- Prohibited for use by individuals with pacemakers, osteoporosis, loss of temperature or pain sensitivity, or skin diseases.
- Avoid physical impacts and dropping the device.

Disposal

Do not dispose of this product with household waste. Please take it to a specialized recycling center to help protect the environment.

Warranty

The MEI1 Face and neck EMS massager by WAU is covered by a 1-year warranty from the date of purchase. If technical defects are found within the warranty period, the device will be replaced or repaired. This warranty does not cover malfunctions caused by improper use or transportation, or devices that have undergone self-repair.

MEI1

FACE AND NECK EMS MASSAGER

Model: **BPWE-881**

Material: **ABS plastic**

Dimensions (HxWxD): **140x80x10 mm**

Weight: **78 g**

Voltage: **5 V**

Power: **4,4 W**

Origin: **PRC**

Authorized manufacturer representative: WAU Global Limited
63A, College road, Northcote, Auckland, 0627, New Zealand
e-mail: info@wauglobal.com

Important Notice: to ensure continuous product improvement,
the manufacturer reserves the right to modify design
and technical specifications. The latest version of this manual
is available at: wauglobal.com

RoHS